



IHE Integration Statement			
Vendor	Product Name	Version	Date
ScImage, Inc.	PicomEnterprise	3.3.x	August 1, 2020
This product implements all transactions required in the IHE Technical Framework to support the IHE Integration Profiles, Actors and Options listed below:			

Cardiology Profiles		
Profiles Implemented	Actors Implemented	Options Implemented
CATH Cardiac Catheterization Workflow	Image Manager/Image Archive Image Display Performed Procedure Step Manager Evidence Creator	Cardiac Cath, Intermittently Connected Modality
ECHO Echocardiography Workflow	Image Manager/Image Archive Image Display Performed Procedure Step Manager Evidence Creator	Echocardiography Intermittently Connected Modality Stress Echo
ECG Retrieve ECG for Display	Image Display Information Source	
ED Evidence Documents	Image Manager/Image Archive Image Display Report Creator	

IT Profiles		
Profiles Implemented	Actors Implemented	Options Implemented
ATNA Audit Trail and Node Authentication	Secure Application Secure Node	
CT Consistent Time	Time Client	SNTP
PDQ Patient Demographics Query	Patient Demographics Consumer	
PIX Patient Identifier Cross Referencing	Patient Identifier Cross Reference Manager	
RID Retrieve Information for Display	Display	
XDS.b Cross Enterprise Document Sharing	Document Consumer	
XDS-SD Cross Enterprise Sharing of Scanned Documents	Document Consumer	

Radiology Profiles		
Profiles Implemented	Actors Implemented	Options Implemented
ARI Access to Radiology Information	Image Manager/Image Archive Image Display Report Reader	Multiple Sources



	Report Repository External Report Repository Access	
CPI Consistent Presentation of Images	Image Manager/Image Archive Image Display Evidence Creator	
DBT Digital Breast Tomosynthesis	Image Manager/Image Archive Image Display	
ED Evidence Documents	Image Manager/Image Archive Image Display Report Creator	
IOCM Image Object Change Management	Image Manager/Image Archive Image Display	
IRWF Import Reconciliation Workflow	Importer Image Manager/Image Archive Performed Procedure Step Manager	
KIN Key Image Note	Image Manager/Image Archive Image Display Evidence Creator	
MAMMO Mammography Image	Image Manager/Image Archive Image Display Evidence Creator	
NM Nuclear Medicine Image	Image Manager/Image Archive Image Display	
PDI Portable Data for Imaging Integration	Portable Media Creator Portable Media Importer Image Display Report Reader	
PGP Presentation of Grouped Procedures	Image Manager/Image Archive Performed Procedure Step Manager	
PIR Patient Information Reconciliation	Image Manager/Image Archive Performed Procedure Step Manager	
PWF Post-Processing Workflow	Image Manager/Image Archive Image Display	
RWF Reporting Workflow	Image Manager/Image Archive Performed Procedure Step Manager Report Manager	
SINR Simple Image and Numeric Report	Report Creator Report Manager Report Repository Report Reader	
SWF Scheduled Workflow	Image Manager/Image Archive Performed Procedure Step Manager Image Display Evidence Creator	
TCE Teaching File and Clinical Trial Export	Export Selector Export Manager	
XDS-I.b Cross-Enterprise Document Sharing for Imaging	Imaging Document Source Imaging Document Consumer	PDF Report Retrieve Images Retrieve Evidence Documents

ScImage's IHE Information and Related Standards Conformance Statements		
HL7	scimage.com/HL7	
DICOM	scimage.com/DICOM	
General IHE Information		
North America: www.ihe.net	Europe: www.ihe-europe.net	Japan: www.ihe-j.org